Instructions for use .....	2
Mode d'emploi.....	5
Istruzioni per l'uso .....	8
Gebrauchsanleitung .....	11
Instrucciones de uso .....	14
Brugsanvisning .....	17
Οδηγίες χρήσης .....	20
Gebruiksaanwijzing .....	23
Bruksanvisning .....	26
Bruksanvisning .....	29
Kullanma Talimatları.....	32
Uputstvo za upotrebu .....	35

## INTENDED USE/INDICATIONS FOR USE

The LINX® Reflux Management System is intended for use in the treatment of symptoms associated with Gastroesophageal Reflux Disease (GERD). The LINX® device is placed at the area of the Lower Esophageal Sphincter (LES) designed to augment a weak LES and minimize or eliminate GERD related symptoms.

The Torax LINX® device is indicated for those subjects diagnosed with pathologic Gastroesophageal Reflux Disease (GERD) defined by abnormal pH testing and who continue to have chronic GERD symptoms despite maximum medical therapy.

The LINX® Reflux Management System is labeled for use by physicians only.

## SYSTEM DESCRIPTION

The LINX® Reflux Management System consists of a series of titanium beads with magnetic cores that are connected with independent titanium wires to form an annular shape, reference Figures 1 and 2.


Figure 1 – Implant at Area of LES


Figure 2 - Illustration of Implant

The LINX® Reflux Management System is comprised of the following components:

- LINX® Device
- LINX® Esophagus Sizing Tool

## CONTRAINDICATIONS

1. Do not implant the LINX® Reflux Management System in patients with suspected or known allergies to titanium, stainless steel, nickel, or ferrous materials. Please note nickel is not a metal component of the LINX® implant, only the sizing tool which is a single-use instrument used only during surgical placement of the implant.

## POTENTIAL COMPLICATIONS

The following is a list of potential complications that may occur with the implantation of the LINX® device: These may include, but may not be limited to the following: Achalasia, Bleeding, Death, Device erosion, Device explant/re-operation, Device failure, Device migration (device does not appear to be at implant site), Diarrhea, Dysphagia, Inability to belch or vomit, Infection, Impaired gastric motility, Injury to the esophagus, spleen, or stomach, Nausea, Odynophagia, Organ damage caused by device migration, Pain, Peritonitis, Pneumothorax, Regurgitation, Saliva/mucus build-up, Stomach Bloating, Vomiting, Worsening of preoperative symptoms (including but not limited to dysphagia or heartburn).

Device erosion refers to at least a portion of the LINX® device passing through the esophageal wall. An observed event has been reported in clinical literature<sup>1</sup>. Patient consultation with a surgeon who is experienced with LINX® procedures and best practices is advised for proper treatment.

## PRECAUTIONS

1. Implantation of the device should only be performed by those physicians who have experience in laparoscopic anti-reflux procedures and has received product specific training.
2. It is the responsibility of the surgeon to advise the patient of the known risks and complications associated with the surgical procedure and implant.
3. The sterile package and device should be inspected prior to use. If sterility or performance of the LINX® device is suspect or compromised, it should not be used.
4. The device is intended for single use only. DO NOT resterilize the device. Functionality and sterility of the device cannot be assured if re-used.
5. The device is magnetic and will be attracted to ferrous objects in the surgical field and other surgical instruments that are ferromagnetic.
6. The LINX® device has not been evaluated in patients with a hiatal hernia >3 cm. Use of the LINX® device in patients with a hiatal hernia >3 cm should be considered on the basis of each patient's medical history and severity of symptoms.

<sup>1</sup> Lipham JC, Taiganides PA, Louie BE, Ganz RA, DeMeester TR. Safety Analysis of first 1000 patients treated with magnetic sphincter augmentation for gastroesophageal reflux disease. *Diseases of the Esophagus* (2014). DOI: 10.1111/dote.12199.

7. Patients should be advised that the LINX® Reflux Management System is a long-term implant. Explant (removal) and replacement surgery may be indicated at any time. Medical management of adverse reactions may include explantation and/or replacement.
8. The safety and effectiveness of the LINX® device has not been evaluated in patients with Barrett's esophagus or Grade C or D (LA classification) or Grade IV (Savary-Miller) esophagitis.
9. The safety and effectiveness of the LINX® device has not been evaluated in patients with electrical implants such as pacemakers and defibrillators, or other metallic, abdominal implants.
10. The safety and effectiveness of the LINX® device has not been evaluated in patients with major motility disorders.
11. The safety and effectiveness of the LINX® Reflux Management System has not been established for the following conditions:
  - Scleroderma
  - Suspected or confirmed esophageal or gastric cancer
  - Prior esophageal or gastric surgery or endoscopic intervention
  - Distal esophageal motility less than 35 mmHg peristaltic amplitude on wet swallows or <70% (propulsive) peristaltic sequences or a known motility disorder such as Achalasia, Nutcracker Esophagus, and Diffuse Esophageal Spasm or Hypertensive LES
  - Symptoms of dysphagia more than once per week within the last 3 months
  - Esophageal stricture or gross esophageal anatomic abnormalities (Schatzki's ring, obstructive lesions, etc.)
  - Esophageal or gastric varices
  - Lactating, pregnant or plan to become pregnant
  - Morbid obesity (BMI >35)

## WARNINGS

1. The device is to be placed around the esophagus including the anterior and excluding the posterior vagus nerve bundle. The device should never be placed outside both vagus nerve bundles.
2. Failure to secure the LINX® device properly may result in its subsequent displacement and necessitate a second operation.
3. The device contains magnets. The LINX® device is considered MR Conditional. Please refer to the MR Conditional Guidelines below for additional detail. Exposure to an MRI environment above the MR Conditional guidelines could cause serious injury to the patient and/or interfere with the magnetic strength and the function of the device. The device is NOT safe when exposed to MRI magnetic fields greater than 0.7 T. In the event alternative diagnostic procedures cannot be used and MRI is required beyond the MR Conditional guidelines, the LINX® device can be safely removed utilizing a laparoscopic technique that does not compromise the option for traditional anti-reflux procedures.
4. The packaged device should not be exposed to temperatures above 60 degrees C as this could adversely affect the magnets and the function of the package/device.

## MRI SAFETY INFORMATION

**MR Conditional**

Non-clinical testing has demonstrated that the LINX® System is MR Conditional. This device can be scanned safely under the following conditions:

- Static magnetic field 0.7-Tesla (0.7 T)
- Spatial gradient field of up to 364 G/cm (3.64 T/m)
- Maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg for 15 minutes of scanning in First Level Controlled Mode
- The LINX® System contains permanent magnets. Should the patient experience discomfort or pain, immediately discontinue the scan and remove the patient from the MR environment.

## RF heating

In non-clinical testing, with body coil excitation, the LINX® System produced a temperature rise of 0.6°C at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 0.7 T Hitachi Altair Open MR scanner with software version 5.1H-1.

Caution: The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

## MR Artifact

In testing using a 0.7 T system with gradient-echo sequencing, the shape of the image artifact follows the approximate contour of the device and extends radially up to 10.3 cm from the implant.

## DIRECTIONS FOR USE

1. Surgical Access:
  - 1.1 Gain surgical access through a laparoscopic port to the esophagus at the region of the gastroesophageal junction.
  - 1.2 Dissect the soft tissues away from the outside of the esophagus at the location of the gastroesophageal junction. Tissue should be removed to expose the outer muscle of the esophagus. Create a tunnel under the posterior vagus nerve through the peri-neural tissue. The anterior vagus nerve will be included within the implant. Care should be taken to avoid injuring the vagus nerve bundles.
2. Sizing of the Esophagus:
  - 2.1 Bring sizing tool into surgical field through a laparoscopic port of appropriate internal diameter.
  - 2.2 Place the sizing tool around the esophagus in the dissected space around the exposed outer muscle and through the tunnel created under the posterior vagus nerve bundle.
  - 2.3 Perform sizing per the appropriate sizing tool instructions for use.
3. Placement of the Implant:
  - 3.1 Bring the chosen device into the surgical field through a laparoscopic port of minimum internal diameter of 10 mm.
  - 3.2 Place the device around the esophagus in the same location that was measured, reference Figure 3.
  - 3.3 Bring the device ends together, the ends are magnetic and will attract when brought in proximity.
  - 3.4 While grasping at the distal knots align and mate the interlocking ends, reference Figure 4.


Figure 3 – Implant Around Esophagus


Figure 4 – Ends Interlock

- 3.5 Make sure the ends are fully engaged. It may be necessary to fully close the ends using a grasper, reference Figure 5.
- 3.6 Trim the suture ends below the proximal knot and remove all excess suture material, reference Figure 6.


Figure 5 – Ends Fully Engaged


Figure 6 – Trimmed Suture Ends

## PACKAGING/STORAGE

The LINX® device is provided sterile and designed to remain sterile unless the primary product pouch has been opened or damaged. Store in a cool, dry place. If opened and not used, discard device or return device to Torax Medical Inc. Do Not Resterilize.

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